

Viral Nucleic Acid Extraction Protocol

Important Notes:

- Wear gloves and lab coat when handling the buffers provided in this kit. They contain irritants.
 - Do not add acidic solution or bleach directly to the sample preparation waste.
 - The precipitates that appear during storage, especially at low temperature, in some buffers are easily dissolved by incubating the bottles at 60°C.
 - Store Proteinase K at -20°C.
 - Add 20ml absolute ethanol to Wash Buffer I and 40ml to Wash Buffer II.
 - * Set thermoblock or water bath at 65°C.
 - * Preheat Elution Buffer at 65°C and at the end, use this warm buffer.
 - * Make sure that ethanol has been added into wash buffers.
- 1- Add 200µl **Binding Buffer** and 10µl **proteinase K** to the microtube containing 200µl of sample. Vortex and spin briefly and incubate at 65°C for 15 min.
 - 2- Add 100µl isopropanol then vortex and spin briefly. Add 1-10µl nucleic acid carrier if needed.
 - 3- Transfer the mixture into the spin column placed in a collection tube.
 - 4- Centrifuge the sample for 1 min at 10000 rpm.
 - 5- Discard flow-through and place the column in a new collection tube.
 - 6- Add 500µl **Wash Buffer I** to the column and centrifuge for 1 min at 10000 rpm.
 - 7- Discard flow-through. Place the column in a new collection tube, add 500µl **Wash Buffer II**, centrifuge for 1 min at 10000 rpm. Then proceed centrifugation at full speed (14000 rpm) for 2-3 min to dry the membrane completely.
 - 8- Place the column in a clean 1.5ml microtube. Open the lid of the column and apply 50-100µl preheated **Elution Buffer** to the center of the membrane.
 - 9- Close the lid and incubate at 65°C for 2-5 min. Centrifuge at full speed (>12000 rpm) for 1 min.
 - 10- Store eluted nucleic acid at -20°C or -70°C.