

Viral Nucleic Acid Extraction

1 Sample (200 μ l)
Binding Buffer (200 μ l)
Proteinase K (10 μ l)

Incubate:
65 °C/15 min

Cent:
1 min
10000 rpm

2 Isopropanol (100 μ l)
Nucleic acid carrier (1-10 μ l)
Transfer to spin column

3 Change collection tube
Wash Buffer I (500 μ l)

Cent:
1 min
10000 rpm

Cent:
1 min/10000 rpm
2-3min/14000 rpm

4 Change collection tube
Wash Buffer II (500 μ l)

5 Place spin column in 1.5ml microtube
Elution Buffer (50-100 μ l)

Incubate:
65 °C/2-5 min
Cent:
1min/12000