

NamaZol RNA Extraction

1 NamaZol (750 μ l)
Tissue (5-50mg) / Blood (250 μ l)

Incubate:
RT
5-15 min

Incubate:
RT/2-5 min
Cen:
15 min/12000 rpm/ 4 °C

2 Chloroform (200 μ l)
Shake for 15 sec

3 Transfer upper phase to new tube
Isopropanol (500 μ l)

Incubate:
RT/5-10 min
Cen:
10 min/12000 rpm/ 4 °C

Cen:
5min/8000 rpm/4 °C

4 Ethanol 75% (1ml)

5 Remove Ethanol & air-drying
DEPC-trated water (20-100 μ l)

Incubate:
55 °C
10-15 min